

EQUINOX AESTHETIC ACADEMY

**OFFICIAL AUSTRALIAN SUPPLIER OF
LIFSKIN PDO THREADS**

STERILE SINGLE USE ABSORBABLE
POLYDIOXANONE SUTURES WITH NEEDLE





AUSTRALIAN SPONSOR

29 Harmony Avenue

Atwell, WA 6164 Australia

Phone: +61 426 969 719

Website: thread.equinox.aestheticacademy.com.au

Email: training@threadlift.com.au

MANUFACTURER

- MECOBI Co Ltd
- Address: 3-101 102 107 108 42-10
Taejanggongdan-gil Wonju-si Gangwon-do
26311 South Korea
- Contact Information:
- Phone: +82 33 733 1612
- Fax: +82 33 733 1614
- Website: www.ducosmedi.com

CE 2195 ARTG KFDA




DOUBLE SCREW



TWO TWISTED SCREW THREAD

Needle		Thread		Model Code
Gauge	Length	USP	Length	
30G	25mm	6-0	30mm	SC30025710NS
	38mm		50mm	SC30038710NS
29G	25mm	6-0	30mm	SC29025610NS
	38mm		50mm	SC29038610NS
	50mm		70mm	SC29050610NS
	60mm		90mm	SC29060610NS
27G	38mm	5-0	50mm	SC27038510NS
	50mm		70mm	SC27050510NS
	60mm		90mm	SC27060510NS
26G	38mm	5-0	50mm	SC26038510NS
	50mm		70mm	SC26050510NS
	60mm		90mm	SC26050510NS
	90mm		150mm	SC26090510NS

- Maximising the effect of screw type. It gives a more powerful effect than a single screw.
- Packing Unit : 20ea/pack, 5packs/box (100ea).



CE 2195 ARTG

SCREW




SINGLE THREAD

Needle		Thread		Model Code
Gauge	Length	USP	Length	
30G	25mm	6-0	30mm	SC30025710NS
	38mm		50mm	SC30038710NS
29G	25mm	6-0	30mm	SC29025610NS
	38mm		50mm	SC29038610NS
	50mm		70mm	SC29050610NS
	60mm		90mm	SC29060610NS
27G	38mm	5-0	50mm	SC27038510NS
	50mm		70mm	SC27050510NS
	60mm		90mm	SC27060510NS
26G	38mm	5-0	50mm	SC26038510NS
	50mm		70mm	SC26050510NS
	60mm		90mm	SC26050510NS
	90mm		150mm	SC26090510NS

- Maximizing the effect of screw type. It gives more powerful effect than single mono.

Screw (Smooth thread twists around the needle, which gives a stronger stimulus to tissues than a straight formed thread), and Double screw (the same as Screw type, but with two smooth threads spirally loaded together).

- Packing Unit: 20 ea/pack, 5 packs/box (100 ea).



CE 2195 ARTG

TWIN THREAD



TWO TWISTED THREAD

Needle		Thread		Model Code
Gauge	Length	USP	Length	
27G	38mm	7-0 x 2	50mm	SC27038710W
	50mm		70mm	SC27050710W
	60mm		90mm	SC27060710W
26G	38mm	7-0, 6-0	50mm	SC26038610W
	50mm		70mm	SC26050610W
	60mm		90mm	SC26060610W
	90mm		150mm	SC26090610W
25G	60mm	6-0 x 2	90mm	SC25060510W
	90mm		150mm	SC25090510W

- Twisted threads are wound around the needle.
- Maximising the effect of screw type. It gives a more powerful effect than a single screw.
- Highly recommended for lifting, tightening, and wrinkle removal.
- CE, ISO GMP approved.
- Packing Unit: 20 ea/pack, 5 packs/box (100 ea).



LIFSKIN



KFDA

CE 2195 ARTG

TRIPLE THREAD



THREE TWISTED THREAD

Needle		Thread		Model Code
Gauge	Length	USP	Length	
26G	38mm	7-0 x 3	50mm	SC26038710R
	50mm		70mm	SC26050710R
	60mm		90mm	SC26060710R
	90mm		150mm	SC26090710R
25G	60mm	6-0 x 3	90mm	SC25060610R
	90mm		150mm	SC25090610R
23G	100mm	4-0 x 3	170mm	SC23100410R

- It shows effects in filling up the cheek, chin, and some parts around the eyes, as well as removing wrinkles. It provides increased volume while lifting and tightening the skin over time.
- CE, ISO GMP approved.
- Packing Unit: 20 ea/pack, 5 packs/box (100 ea).



LIFSKIN



KFDA

CE 2195 ARTG



EYE CANNULA

MONO

Needle		Thread		Model Code
Gauge	Length	USP	Length	
29G	38mm	6-0	50mm	SC29038610NE

EYE CANNULA

SCREW

Needle		Thread		Model Code
Gauge	Length	USP	Length	
29G	38mm	6-0	50mm	SC29038610NS

- Mono thread with blunt needle.
- It has been designed to treat eye bags and fine wrinkles around the eyes.
- A blunt needle minimises swelling, bruising, and bleeding during the procedure.

- Packing Unit: 20 ea/pack, 5 packs/box (100 ea).



LIFSKIN



KFDA

CE 2195 ARTG

COG THREAD




SHARP NEEDLE TYPE (SEMI - BLUNT NEEDLE)

4D, 360°, BI-DIRECTION

Needle		Thread		Model Code
Gauge	Length	USP	Length	
19G	70mm	1-0	90mm	SCD190700U
	100mm		150mm	SCA191000C
21G	60mm	2-0	90mm	SCA210620C
	90mm		150mm	SCA210920C
23G	60mm	3-0	90mm	SCA230630C
	90mm		150mm	SCA230930C

360°, bi-directional cog thread with sharp needle, 45° blunt, L type and You type blunt needle • Mostly used for face, eyebrows and body lifting • CE, ISO, GMP approved

- Packing Unit: 10 ea/pack, 10 packs/box (100 ea).



LIFSKIN



KFDA

CE 2195 ARTG

COG THREAD




360°, BI-DIRECTION

L-TYPE CANNULA

Needle		Thread		Model Code
Gauge	Length	USP	Length	
18G	100mm	2	150mm	SCB181000U
19G	70mm	1-0	90mm	SCB190700U
	100mm		150mm	SCB191000C
21G	60mm	2-0	90mm	SCB210620C
	90mm		150mm	SCB210920C
23G	60mm	3-0	90mm	SCB230630C
	90mm		150mm	SCB230930C

360°, bi-directional Cog thread with sharp needle, 45° blunt, L type and U type blunt needle

- Mostly used for face, eyebrows and body lifting
- CE, ISO, GMP approved
- Packing Unit: 10 ea/pack, 10 packs/box (100 ea).



LIFSKIN



KFDA

CE 2195 ARTG

COG THREAD




U-TYPE CANNULA

360°, BI-DIRECTION

Needle		Thread		Model Code
Gauge	Length	USP	Length	
18G	100mm	2	150mm	SCB181000U
19G	70mm	1-0	90mm	SCB190700U
	100mm		150mm	SCB191000C
21G	60mm	2-0	90mm	SCB210620C
	90mm		150mm	SCB210920C
23G	60mm	3-0	90mm	SCB230630C
	90mm		150mm	SCB230930C

- Bi-directional 4D cutting Cog thread with U-type blunt needle.
- Specially designed to strongly lift sagging skin and improve wrinkles, applicable on the face and body.
- The best item of Cog threads.
- Packing Unit: 10 ea/pack, 10 packs/box (100 ea).



LIFSKIN



KFDA

CE 2195 ARTG

NOSE




L-TYPE CANNULA

Needle		Thread		Model Code
Gauge	Length	USP	Length	
19G	38mm	2	50mm	SCB190380C
	50mm		80mm	SCB190500C

- Cog thread with L-type cannula needle. • It has been designed to treat the nose without incision, providing an easy and safe treatment!
- Packing Unit: 10 ea/pack, 10 packs/box (100 ea).

The nose thread series is specially designed for nose procedures of stereoscopic volume and holding.



LIFSKIN



KFDA

CE 2195 ARTG

MOLDING COG




W-TYPE CANNULA

BI-DIRECTION

Needle		Thread		Model Code
Gauge	Length	USP	Length	
18G	100mm	2	150mm	SCD191000C
19G	100mm	3	150mm	SCC181000U

- Bi-direction 4D cutting Cog thread with W type blunt needle. • Specially designed to strongly lift sagging skin and improve wrinkles, Applicable on face and Body. • **The best item of Cog threads.** •
- Packing Unit : 10ea/pack, 10packs/box(100ea).

The advantage of this product is the strength and durability of the thread. Unlike other types of cog thread products, this type uses a manufacturing method called the “moulding press process,” which results in very minor force loss in comparison with other types that are produced by cog pattern processing.



LIFSKIN



KFDA

CE 2195 ARTG

Multi Thread



L-TYPE CANNULA

Needle		Thread		Model Code
Gauge	Length	USP	Length	
19G	38mm	2	50mm	SCBI90380C
	50mm		80mm	SCBI90500C

. 10 Mono threads in one needle. Powerful volume-up effect. Much more collagen is collected inside the L-type cannula needle: It minimises swelling, bruising, and bleeding during the procedure.

- Packing Unit: 10 ea/pack, 10 packs/box (100 ea)

The Multi-Fill product is made up of many strands with a blunt cannula. Similar to the Cavern-Fill, it interacts with Cavern-Fill and is suitable for use in the areas of the forehead, chin, lateral cheek hollow, lower lid, tear trough, marionette line, and lip augmentation. Additionally, it can be used for eyebrow augmentation, as well as for nose and neck wrinkles. It has a mesh-like shape so that connective tissue composed of collagen is effectively formed around the thread. Unlike the Cog thread, it acts as a volume supplement rather than having a physical pulling effect.



LIFSKIN



KFDA

CE 2195 ARTG


Double Needle Combi



360, BI-DIRECTION

Needle		Thread		Model Code
Gauge	Length	USP	Length	
19G	110mm	1-0	290mm	SCC191100US
			430mm	SCC191100UL

This high-end product consists of one suture (300mm or 420mm) and two needles which are specially designed for professional use. Dual needle products have their own advantages, in that users can manually control with accuracy the detailed procedure pathway, which can also result in a symmetrical performance. Professionals recommend blunt tip needles to reduce damage to the blood vessels and to minimise bruising during insertion. Dual needle threads are typically used in the forehead, lateral face, medial area, mouth corner, etc.



LIFSKIN



KFDA

CE 2195 ARTG

EQUINOX AESTHETIC ACADEMY

**OFFICIAL AUSTRALIAN SUPPLIER OF
LIFSKIN PDO THREADS**



Equinox training Services Pty Ltd

ACN: 643865543

29 Harmony Avenue

Atwell, WA 6164 Australia

Phone: +61 426 969 719

Website: thread.equinoxestheticacademy.com.au

Email: training@threadlift.com.au